




MagCapture[™] EXOSOME ISOLATION KIT PS Ver.2

EXOSOME ISOLATION BY A NOVEL AND IMPROVED AFFINITY METHOD

INTRODUCTION

MagCapture[™] Exosome Isolation Kit PS Ver.2 adopts a novel affinity purification method using magnetic beads and a phosphatidylserine (PS)-binding protein. By using a phosphatidylserine (PS)-binding protein, extracellular vesicles are captured in a metal ion-dependent manner and are subsequently eluted from magnetic beads with a metal-chelating reagent at neutral pH.



Features

- + Isolation by PS affinity method
- + Ultracentrifugation is not required
- + High purity and high yield of intact EVs
- + Highly reproducible yield
- + Improved recovery rate
- + Short operation time (~1.5 hours)
- + No preservatives (less cytotoxicity)

Applicable Samples

- + Cell culture supernatant
- + Serum, plasma
- + Urine, saliva, etc

Kit Component

		2 TESTS	10 TESTS
1	Biotin Capture Magnetic Beads	120 μ L \times 1 tube	$600 \mu L \times 1 \text{tube}$
2	Biotin-labeled Exosome Capture	20 μ L \times 1 tube	100 μL × 1 tube
3	Exosome Immobilizing/Washing Buffer (10x)	5 mL × 1 bottle	25 mL × 1 bottle
4	Exosome Binding Enhancer (500×)	300 μL \times 1 tube	1500 μL × 1 tube
5	Exosome Elution Buffer (10x)	$300~\mu L \times 1~tube$	1500 μ L \times 1 tube
6	Reaction Tubes	4 tubes	22 tubes

ORDERING INFORMATION

WAKO CODE	PRODUCT NAME	STORAGE CONDITION	PACKAGE SIZE
294-84101	MagCapture™ Exosome Isolation Kit PS Ver.2 NEW	Keep at 2-10 °C	2 Tests
290-84103	MagCapture™ Exosome Isolation Kit PS Ver.2 NEW	Keep at 2-10 °C	10 Tests

COMPARISON WITH OTHER PURIFICATION METHODS


	METHOD	EXOSOME PURITY	EXOSOME RECOVERY	INTACT VESICLES RECOVERY
	PS affinity method	++++	++++	Yes
I	Ultracentrifugation	+ +	+ +	Yes
	Polymer-based precipitation	+	+ +	Yes
	Exosome surface antigen affinity method (using antibodies)	++++	++	No
	+ Very low + + Low + + + Medium + + + + + High + + + + + Very High			


COMPARISON WITH CURRENT EV ISOLATION KIT

Γ		CURRENT KIT	VER. 2
	Recovery amount (culture supematant & blood sample	+++	++++
	Reaction time	3 hour +	1 hour +
	Reusability (Magnetic beads)	5 times	5 times
	Cell cytotoxicity (Elution buffer)	depending on cell lines	low
	+ Very low + + Low + + + Medium + + + +	High +++++ Very High	


EVs from MSC culture supernatant were isolated by our current kit (code: 299-77603, 293-77601) and the new Ver.2 kit and examined by TEM, nanoparticle tracking analysis (NTA), and Western blot.

(1) EV morphology







EVs with more round and uniform morphology could be obtained via isolation by the new Ver.2 kit.



(3) Cell cytotoxicity

EVs from COLO201 cell culture supernatant were isolated and purified by the current kit and Ver.2 kit. Same amount of purified EVs and the exosome elution buffer were added into human normal fibroblasts cells.





After 48 hours, cell morphological changes and cell death were observed with the current kit. In contrast, these effects could not be obtained for Ver.2 kit. Ver.2 kit is applicable for uptake experiments (in vitro or in vivo) without buffer exchange.

REFERENCE

1. A novel affinity-based method for the isolation of highly purified extracellular vesicles, W. Nakai, T. Yoshida, D. Diez, Y. Miyatake, T. Nishibu, N. Imawaka, K. Naruse, Y. Sadamura and R. Hamayama, Sci. Rep., 6, 33935 (2016).

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit our online catalog to search for other products from FUJIFILM Wako; https://labchem-wako.fujifilm.com / This leaflet may contain products that can not be exported to your country due to regulations. / Bulk quote requests are welcome. Please contact us.

DISTRIBUTION & INQUIRIES

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstr. 12 41468 Neuss · Germany

Phone +49 2131 311 271 Fax +49 2131 311 110

labchem_wkeu@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road Richmond VA 23237 · USA

Phone +1-804-714-1920 Fax +1-804-271-7791 wkuslabchem@fujifilm.com

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome Chuo-Ku Osaka 540-8605 · Japan

Phone +816 6203 3741 Fax +81 6 6203 1999 ffwk-cservice@fujifilm.com

